

Software Design Notes

1. Revision History

Rev.	Revision Notes	DCN	Date
1	Initial Release.	N/A	24SEP2021
2			

2. Introduction

2.1. Purpose

This document presents software design notes for VBA code that logs a form and saves the form as a PDF.

This software is intended to automatically log and file the Ultrasound Transducer Test Form, FRM.100-003 when the form is printed.

2.2. Scope

This document applies to the VBA subroutine(s) listed in Table 1 – VBA Form Submission Subroutine.

Table 1 – VBA Form Submission Subroutine

Subroutine	Description
Document_Open()	<ul style="list-style-type: none">Runs when the form document is opened.Located in ThisDocument.cls
Register_Event_Handler()	<ul style="list-style-type: none">Runs when called.Located in Reg_Event_Handler.basCreates an instance of the Word Application object so that application events can be used.
wordApp_DocumentBeforePrint()	<ul style="list-style-type: none">Runs when DocumentBeforePrint application event is triggered.Located in EventClassModule.clsLogs and creates PDF of form.

2.3. References

2.3.1. Files

File	Location
Form	C:\Repos\Software\Bexar VBA\BexarVBA\100\003 Ultrasound Transducer Test Form
Log	C:\Repos\Software\Bexar VBA\BexarVBA\100\003 Ultrasound Transducer Test Form\data
PDFs	C:\Repos\Software\Bexar VBA\BexarVBA\100\003 Ultrasound Transducer Test Form\pdfs
VBA	https://github.com/nolan210/BexarVBA.git

3. Software Specification Requirements

SSR1 – Usability

No text (title)

SSR1.1

Published in a Word Macro Enabled Document.

SSR1.2

Accesses and saves files on a local file server.

SSR1.3

Does not require the document user to save the Word document.

SSR1.4

All logging and PDF saving functionality is accomplished automatically, without user input, before the document is printed.

SSR2 – Functionality

No text (title)

SSR2.1

VBA Project (title)

SSR2.1.1

Protect the VBA source code with the password in Table 2 – VBA Protection Password.

Table 2 – VBA Protection Password

Password	9hvmq6KchrZ+V&&j
-----------------	------------------

SSR2.2

Log Data (title)

SSR2.2.1

Log the date and time at print event in the Excel file.

SSR2.2.2

Log the username of account logged into computer at print event in the Excel file.

SSR2.2.3

Log the computer ID at print event in the Excel file.

SSR2.2.4

Log the form data in the Excel file.

SSR2.3

Excel File (title)

SSR2.3.1

Excel file with a log sheet at the file location in Table 3 – Excel File Location.

Table 3 – Excel File Location

File Location	C:\Repos\Software\Bexar VBA\BexarVBA\100\003 Ultrasound Transducer Test Form\data\Log.xlsx
----------------------	--

SSR2.3.2

Protect the Excel file's data sheet with the password in Table 4 – Excel Sheet Protection Password.

Table 4 – Excel Sheet Protection Password

Password	WT\$?@L&x=F5qYfC#
-----------------	-------------------

SSR2.4

PDFs (title)

SSR2.4.1

Save PDF of the form being printed according to details in Table 5 – Print PDF Details.

Table 5 – Print PDF Details

File Location	C:\Repos\Software\Bexar VBA\BexarVBA\100\003 Ultrasound Transducer Test Form\data\pdfs\
File Name	YYYYMMDDHHMMSS_username.pdf

4. Source Code

Table 6 – ThisDocument.cls VBA Source Code

```
VERSION 1.0 CLASS
BEGIN
    MultiUse = -1    'True
END
Attribute VB_Name = "ThisDocument"
Attribute VB_GlobalNameSpace = False
Attribute VB_Creatable = False
Attribute VB_PredeclaredId = True
Attribute VB_Exposed = True
Private Sub Document_Open()
    Register_Event_Handler
End Sub
```

Table 7 – Reg_Event_Handler.bas VBA Source Code

```
Attribute VB_Name = "Reg_Event_Handler"
Dim x As New EventClassModule
Sub Register_Event_Handler()
    Set x.wordApp = Word.Application
End Sub
```

Table 8 – EventClassModule.cls VBA Source Code

```
VERSION 1.0 CLASS
BEGIN
    MultiUse = -1    'True
END
Attribute VB_Name = "EventClassModule"
Attribute VB_GlobalNameSpace = False
Attribute VB_Creatable = False
Attribute VB_PredeclaredId = False
Attribute VB_Exposed = False
'Nolan Manteufel
'Document No. 860-113-980
'Document Rev. 00
'15SEP2020
'VBA Password = 9hvmq6KchrZ+V&&j

'Used to move documents
Private Const PDF_FILEPATH As String = "O:\data\INF970-130\"
Private Const LOG_FILEPATHNAME As String = "O:\data\INF970-130\Log.xlsx"
Private Const LOG_PASSWORD As String = "Z45KUWkf37Drdy*2"

'Used to access CC form fields
Private Const ccSN = "ccSN"
Private Const ccCondition = "ccCondition"
Private Const ccTester = "ccTester"
Private Const ccUC = "ccUC"
Private Const ccTicket = "ccTicket"
Private Const ccIssue = "ccIssue"
Private Const ccT1OpMode = "ccT1OpMode"
Private Const ccT1Param = "ccT1Param"
Private Const ccT1Start = "ccT1Start"
Private Const ccT1Stop = "ccT1Stop"
Private Const ccT1Marker = "ccT1Marker"
Private Const ccT2OpMode = "ccT2OpMode"
Private Const ccT2Frequency = "ccT2Frequency"
Private Const ccYmax = "ccYmax"
Private Const ccYmin = "ccYmin"
Private Const ccMarkerReading = "ccMarkerReading"
```

```
Private Const ccVSWR = "ccVSWR"
Private Const ccConclusion = "ccConclusion"

'Used to access CC document stamps
Private Const ccUser = "ccUser"
Private Const ccComputer = "ccComputer"
Private Const ccDatetime = "ccDatetime"

Private Const NUMBER_OF_DATA_FIELDS = 21

'Used to access data array and determine data locations in log
Enum formInputEnums
    iUser
    iComputer
    iDatetime
    iSN
    iCondition
    iTester
    iUC
    iTicket
    iIssue
    iT1OpMode
    iT1Param
    iT1Start
    iT1Stop
    iT1Marker
    iT2OpMode
    iT2Frequency
    iYmax
    iYmin
    iMarkerReading
    iVSWR
    iConclusion
    iPDF
End Enum

Public WithEvents wordApp As Word.Application
Attribute wordApp.VB_VarHelpID = -1
Private Sub wordApp_DocumentBeforePrint(ByVal Doc As Document, Cancel As Boolean)

'Nolan Manteufel
'15SEP2020

'Logs the form in LOG_FILEPATHNAME excel file.
'Protects the LOG_FILEPATHNAME excel file with LOG_PASSWORD.
'Saves PDF of form in PDF_FILEPATH folder.

'Exit if error
On Error GoTo ErrHandler

'Variables
Dim formDoc As Document
Set formDoc = ActiveDocument
Dim cc As ContentControls

Dim excelApp As Excel.Application
Dim logWB As Excel.Workbook
Dim logSheet As Excel.Worksheet

'User
Dim formUser As String
formUser = Environ$("username")
formUser = LCase(formUser)

'Computer
Dim formComputer As String
formComputer = Environ$("computername")
formComputer = LCase(formComputer)

'Datetime
```

```
Dim dateTime As String
Dim timeOnly As String
timeOnly = Format(Time, "hhmmss")
dateTime = Format(Date, "yyyymmdd")
dateTime = dateTime & timeOnly

'Inputs array
Dim formInputs(NUMBER_OF_DATA_FIELDS) As String
Dim dataColumnIndex As Integer

'Populate metadata
formInputs(iUser) = formUser
formInputs(iComputer) = formComputer
formInputs(iDatetime) = dateTime

'Populate data
formInputs(iSN) = formDoc.SelectContentControlsByTag(ccSN) (1).Range.Text
formInputs(iCondition) = formDoc.SelectContentControlsByTag(ccCondition) (1).Range.Text
formInputs(iTester) = formDoc.SelectContentControlsByTag(ccTester) (1).Range.Text
formInputs(iUC) = formDoc.SelectContentControlsByTag(ccUC) (1).Range.Text
formInputs(iTicket) = formDoc.SelectContentControlsByTag(ccTicket) (1).Range.Text
formInputs(iIssue) = formDoc.SelectContentControlsByTag(ccIssue) (1).Range.Text
formInputs(iT1OpMode) = formDoc.SelectContentControlsByTag(ccT1OpMode) (1).Range.Text
formInputs(iT1Param) = formDoc.SelectContentControlsByTag(ccT1Param) (1).Range.Text
formInputs(iT1Start) = formDoc.SelectContentControlsByTag(ccT1Start) (1).Range.Text
formInputs(iT1Stop) = formDoc.SelectContentControlsByTag(ccT1Stop) (1).Range.Text
formInputs(iT1Marker) = formDoc.SelectContentControlsByTag(ccT1Marker) (1).Range.Text
formInputs(iT2OpMode) = formDoc.SelectContentControlsByTag(ccT2OpMode) (1).Range.Text
formInputs(iT2Frequency) = formDoc.SelectContentControlsByTag(ccT2Frequency) (1).Range.Text
formInputs(iYmax) = formDoc.SelectContentControlsByTag(ccYmax) (1).Range.Text
formInputs(iYmin) = formDoc.SelectContentControlsByTag(ccYmin) (1).Range.Text
formInputs(iMarkerReading) = formDoc.SelectContentControlsByTag(ccMarkerReading) (1).Range.Text
formInputs(iVSWR) = formDoc.SelectContentControlsByTag(ccVSWR) (1).Range.Text
'formInputs(iConclusion) left empty
formInputs(iPDF) = PDF_FILEPATH & dateTime & "_" & formUser & ".pdf"

'Request for confirmation
'Dim intResponse As Integer
'intResponse = MsgBox("Are you sure you wanted to print?", vbYesNo)
'If intResponse = vbNo Then Cancel = True

'Launch excel
Set excelApp = New Excel.Application

'The workbook
Set logWB = excelApp.Workbooks.Open(filename:=LOG_FILEPATHNAME)

'The worksheet
Set logSheet = logWB.Worksheets("log")

'Unprotect the log
logSheet.Unprotect (LOG_PASSWORD)

'Log the entry
Dim logPKey As Long
logPKey = logSheet.UsedRange.Rows.Count + 1
logSheet.Cells(logPKey, 1).Value = logPKey

For dataColumnIndex = 0 To NUMBER_OF_DATA_FIELDS
logSheet.Cells(logPKey, dataColumnIndex + 2).Value = formInputs(dataColumnIndex)
Next

'Protect the log
logSheet.Protect (LOG_PASSWORD)

'Close log
logWB.Close savechanges:=True

'Release object references
Set cc = Nothing
```

```
Set logSheet = Nothing
Set logWB = Nothing
Set excelApp = Nothing

'' STAMPS REMOVED
''Set form field stamps
'formDoc.SelectContentControlsByTag(ccUser)(1).LockContents = False
'formDoc.SelectContentControlsByTag(ccUser)(1).Range.Text = formInputs(iUser)
'formDoc.SelectContentControlsByTag(ccComputer)(1).LockContents = False
'formDoc.SelectContentControlsByTag(ccComputer)(1).Range.Text = formInputs(iComputer)
'formDoc.SelectContentControlsByTag(ccDatetime)(1).LockContents = False
'formDoc.SelectContentControlsByTag(ccDatetime)(1).Range.Text = formInputs(iDatetime)

'Save as PDF
formDoc.ExportAsFixedFormat _
    OutputFileName:=formInputs(iPDF), _
    ExportFormat:=wdExportFormatPDF, _
    OpenAfterExport:=False

''Clear form field stamps
'formDoc.SelectContentControlsByTag(ccUser)(1).Range.Text = " "
'formDoc.SelectContentControlsByTag(ccUser)(1).LockContents = True
'formDoc.SelectContentControlsByTag(ccComputer)(1).Range.Text = " "
'formDoc.SelectContentControlsByTag(ccComputer)(1).LockContents = True
'formDoc.SelectContentControlsByTag(ccDatetime)(1).Range.Text = " "
'formDoc.SelectContentControlsByTag(ccDatetime)(1).LockContents = True

''Confirmation message
'MsgBox ("Awesome, " & UCase(formUser) & "! The form has been saved successfully.")

ErrorHandler:

End Sub
```